






The **ART** of **BEAUTY** combined with the
SCIENCE for giving you results
BEYOND expectations

EIKONHA

Eikonha is next generation cross-linked hyaluronic acid dermal filler with lidocaine designed for architectural reshaping of the Face, based on Enhanced Interlinked 3D Matri (EIM) technology which allows easy to inject with natural-looking, and long-lasting effect.

Art

The artistic approach
through Golden Ratio




Science
EIM Technology


EIM Technology

Eikonha is next generation cross-linked hyaluronic acid dermal filler with lidocaine designed for architectural reshaping of the Face, based on Enhanced Interlinked 3D Matri (EIM) technology which allows easy to inject with natural-looking, and long-lasting effect.


1
Cross-
linking




2
Dialysis



3
Sieving



4
Lubrication



Safety Parameters: Eikonha

<0.25 UI/ml

Low Endotoxin
Level

Endotoxin Level

<2 µg/gm

Non-Detectable
BDDE Residual
Content

BDDE

<6 µg/ml

Low Protein
Load

Protein Load

Benefits to Doctor

-  Higher safety & tolerance
-  Higher stability against the natural degradation process
-  Good cohesivity of gel & easy to inject
-  Indication wise specially designed range
-  Allow nutrition, oxygen, & hormones to go through freely, protect regular physiology of the tissue.

Benefits to Patients

-  Completely safe with no reported adverse effects
-  Quick & long-lasting result
-  Quick & long-lasting result
-  Less painful
-  High satisfaction

GLOW

GLOW

GLOW

GLOW


GLOW

GLOW

GLOW OUT SKIN

EIKONI-HA GLOW

Eikonha Glow is Specially Designed for Filling of Facial Deformaties including Superficial Skin Depressions, Periorbital & Peribuccal Fine Lines, and Forehead Lines.



Indications

- Perioral Lines
- Periorbital Lines
- Superficial Facial Lines
- Skin Depressions
- As a Mesotherapy



HA Concentration: 20mg/ml


Lidocaine Concentration: 0.3%

Volumizing Capacity: +++

Volume: 1 ml

Needles & Cannula: 30½ G needle

EIKONI-HA GLOW



Indications

- Perioral Lines
- Periorbital Lines
- Superficial Facial Lines
- Skin Depressions
- As a Mesotherapy



HA Concentration: 20mg/ml

Lidocaine Concentration: 0.3%

Volume: 1ml

Needles & Cannula: 30½ G needle


POUT



POUT OUT LOUD

EIKONHA POUT

Eikonha Pout is Specially Designed for Filling of Lips to Achieve the Fullness and Pouting Effect. Eikonha Pout Can Also be Injected in the Vermilion Border and for Oral Commissures.



Indications

- Pouting
- Lip Fullness
- Vermilion Border
- Oral Commissures



HA Concentration: **20mg/ml**

Lidocaine Concentration: **0.3%**

Volume: **1ml**

Needles & Cannula: **30½ G needle**

GLAM

GLAM

GLAM

GLAM


GLAM

GLAM

**GLAM OUT WITH
ELEGANCE**

EIKONHA GLAM

Eikonha Glam is Specially Indicated for Deep Depressions and Facial Lines. Inject it into the Deep Dermis for the Contours of the Nose, to Correct the Naso-Labial folds, for Moderate Chin Projection, and for Moderate Cheekbones Augmentation.



Indications

- Nasolabial Folds
- Cheek's Fullness
- Chin Augmentation
- Non-Surgical Rhinoplasty



HA Concentration: 20mg/ml

Lidocaine Concentration: 0.3%

Volume: 1ml

Needles & Cannula: 27½ G needle(2)

PROFI

PRO

PRO

PRO

PR


FI

**FILL OUT WITH
CONFIDENCE**



EIKONHA PROFIL

Eikonha Profile is Specially Designed For Volume Creation by Contouring the Facial Tissue. It also Indicated for Rejuvenation of Hand, Nose bridge, Deep Wrinkles and Folds in Chin and Jawlines.



Indications

- Volume creation
- Facial Contours
- Hands
- Chin & Cheek
- Jawline
- Temple
- Large Deep Wrinkles & Folds



HA Concentration: 20mg/ml

Lidocaine Concentration: 0.3%

Volumizing Capacity: ++++

Volume: 1 ml

Needles & Cannula: 3x23*3/4 Needles
3x25G Cannula

EIKONIHA



Concentration	20mg/ml stabilized hyaluronic acid with Lidocaine
Approximate number of gel particles/ml	300000
Particle size (mm)	0.10 - 0.15
Where to inject	Upper part of dermis
Volume of syringe	1ml
Needle size/ quantity	25G 1"/ 6 needles in pack
Duration	6 Months
Packaging information	1 ml X 1 X 1 blisters

Concentration	20mg/ml stabilized hyaluronic acid with Lidocaine
Approximate number of gel particles/ml	10000
Particle size (mm)	0.28-0.5
Where to inject	Deep layer of dermis and/of Surface layer of subcutis
Volume of syringe	1ml
Needle size/ quantity	27G ½"/ 2 needle in pack
Duration	8 - 12 Months
Packaging information	1 ml X 1 X 1 Blister


Concentration	20mg/ml stabilized hyaluronic acid with Lidocaine
Approximate number of gel particles/ml	100000
Particle size (mm)	0.15 - 0.28
Where to inject	Middle part of dermis
Volume of syringe	1ml
Needle size/ quantity	30G ½"/ 2 needles in pack
Duration	8 - 12 Months
Packaging information	1 ml X 1 X 1 blisters

Concentration	20mg/ml stabilizedhyaluronic acid with Lidocaine
Approximate number of gel particles/ml	5000
Particle size (mm)	0.5 - 1.25
Where to inject	Subcutaneous or Supraperiosteal
Volume of syringe	1ml
Needle size/ quantity	25G 1"/ 6 needles in pack
Duration	12 - 18 Months
Packaging information	1 ml X 1 X 3 blisters

Value added service for your better practice end to end solution for your better services to patient



In-clinic branding




Wall Posters



**Injecting
technique**



**Digital marketing
support**



Patient brochures

EIKONHA



10066



AK Pharma, Inc., 9100 S Dadeland Blvd Suite 1500-MBP# 150006 Miami, FL 33156
web:www.eikonha.com | info@eikonha.com